						Availab	le to					Clanual Convenience to
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
Alignate	Activheal - Alignate/ Alignate rope for Cavity	5 x 5cm	£0.56	٧	٧	٧	٧	٧	٧	Suitable for wide range of wound types particularly		ELS139
Alignate	Activheal - Alignate/ Alignate rope for Cavity	10 x 10cm	£1.07	٧	٧	٧	٧	٧	٧	those with moderate to heavy volumes of exudate. These dressings are not suitable for use on dry	Cut sheets to size of the wound to avoid maceration of the surrounding skin Requires a secondary dressing Por cavity wound only Effective for up to 7days Primary contact layer for a wide variety of wounds Indications include skin tears, burns and surgical wounds For use under leg ulcer bandages	ELS140
Alignate	Activheal - Alignate/ Alignate rope for Cavity	10 x 200cm	£2.74	٧	٧	٧	٧	٧	٧	wounds. Cavity wounds should not be over-packed as the fibres swell and can lead to tissue ischaemia. Cut	-	ELS141
Alignate	Activheal - Alignate/ Alignate rope for Cavity	2.5 x 30cm	£1.74	٧	٧	٧	٧	٧	٧	sheets to the size of the wound to avoid maceration and excoriation of the surrounding skin, as any exudate will transfer laterally.	 Requires a secondary dressing 	ELS142
Gelling Fibre	Durafiber	5 x 5	£0.78	٧	٧	٧	٧	٧	٧	Gelling fibre dressings have the ability to absorb and	. Dominos o consedente decesios	ELY330
Gelling Fibre	Durafiber	10 x 10	£1.85	٧	٧	٧	٧	٧	٧	lock away large amounts of exudate; it prevents lateral wicking of exudate and forms a soft gel. Not suitable for use on dry wounds. Do not use in	Requires a secondary dressing	ELY331
Gelling Fibre	Aquacel ribbon	2x45cm		٧		٧	٧		٧	combination with hydrogels	Requires a secondary dressing Requires a secondary dressing Pror cavity wound only Effective for up to 7days Primary contact layer for a wide variety of wounds Indications include skin tears, burns and surgical wounds For use under leg ulcer bandages E E E E E E E	ELY013
Low adherent wound contact layer	Atrauman-1st choice	5x5cm	£0.78	٧	٧	٧	٧	٧	٧		• Effective for up to 7days	EKB039
Low adherent wound contact layer	Atrauman-1st choice	10 x 10	£1.85	٧	٧	٧	٧	٧	٧	Used to prevent trauma to granulating or friable	, , ,	EKB040
Low adherent wound contact layer	Cuticell Contact	2x45cm		٧	٧	٧	٧	٧	٧	wounds.	· · · · · · · · · · · · · · · · · · ·	
Low adherent wound contact layer	Profore WCL	14x20cm	£0.41		٧	٧			٧			ELA127
Absorbent dressing	Absorbent dressing- Zetuvit E- 1st choice Low to moderate amounts of exudate	10 x 10cm	£0.15	٧	٧	٧	٧	٧	٧			ELZ069
Absorbent dressing	Absorbent dressing- Zetuvit E- 1st choice Low to moderate amounts of exudate	10 x 20cm	£0.19	٧	٧	٧	٧	٧	٧	These are dressing that absorb wound exudate Depending on the amount of execute, low, medi or high levels of absorbency may be required. Cre		EJA026
Absorbent dressing	Absorbent dressing- Zetuvit E- 1st choice Low to moderate amounts of exudate	20 x 20cm	£0.30	٧	٧	٧	٧	٧	٧	and ointments can reduce the absorbency of some of the products. Absorbency may be reduced under compression	Cut sheets to size of the wound to avoid maceration of the surrounding skin Requires a secondary dressing Requires a secondary dressing For cavity wound only Effective for up to 7days Primary contact layer for a wide variety of wounds Indications include skin tears, burns and surgical wounds For use under leg ulcer bandages E E E E E E E E E E E E E	EJA028
Absorbent dressing	Absorbent dressing- Zetuvit E- 1st choice Low to moderate amounts of exudate	20 x 40cm	£0.84	٧	٧	٧	٧	٧	٧			EME050
Absorbent dressing	Kliniderm Superabsorbent- 2nd choice Moderate amounts of execute	10 x 10cm	£0.39		٧	٧			٧			EKH072
Absorbent dressing	Kliniderm Superabsorbent- 2nd choice Moderate amounts of execute	10 x 20cm	£0.56		٧	٧			٧	These are dressing that absorb wound exudate. Depending on the amount of execute, low, medium or high levels of absorbency may be required. Creams		EKH073
Absorbent dressing	Kliniderm Superabsorbent- 2nd choice Moderate amounts of execute	20 x 20cm	£1.11		٧	٧			٧	and ointments can reduce the absorbency of some of the products. Absorbency may be reduced under compression		EKH074
Absorbent dressing	Kliniderm Superabsorbent- 2nd choice Moderate amounts of execute	20 x 30cm	£1.92		٧	٧			٧			EKH075
Superabsorbent dressing	Zetuvit Plus- 1st choice	10 x 10cm	£0.56	٧	٧	٧	٧	٧	٧			EME046
Superabsorbent dressing	Zetuvit Plus- 1st choice	10 x 20cm	£0.79	٧	٧	٧	٧	٧	٧	These are dressings that absorb wound exudate.		EME047
Superabsorbent dressing	Zetuvit Plus- 1st choice	15 x 20cm	£0.85	٧	٧	٧	٧	٧	٧	Superabsorbent has the highest level of absorbency and is well suited to with absorbing high amounts of		EME048
Superabsorbent dressing	Zetuvit Plus- 1st choice	20 x 25cm	£1.19	٧	٧	٧	٧	٧	٧	exudate.		EME049
Superabsorbent dressing	Zetuvit Plus- 1st choice	20 x 45cm	£1.92	٧	٧	٧	٧	٧	٧			EME050

D 1 01		C.				Availab	le to			0.77.11.6		
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
Superabsorbent dressing	Eclypse- High exuding wounds	10 x 10cm	£0.77		٧	٧			٧	These are dressings that absorb wound exudate.		EJE029
Superabsorbent dressing	Eclypse- High exuding wounds	15 x 15cm	£1.02		٧	٧			٧	Superabsorbent has the highest level of absorbency and is well suited to with absorbing high amounts of		EJE032
Superabsorbent dressing	Eclypse- High exuding wounds	20 x 30cm	£2.20		٧	٧			٧	exudate.		EJE030
Form dressing	Mepilex Border- 1st choice	7 x 7.5cm	£1.26	٧	٧	٧	٧	٧	٧	Facus dusseines and he would be absorb low to bigh		ELA822
Form dressing	Mepilex Border- 1st choice	10 x 10.5cm	£1.68	٧.	٧	٧	٧	٧	٧	Foam dressings can be used to absorb low to high		ELA828
Form dressing	Mepilex Border- 1st choice	10 x 12.5cm	£2.44	٧	٧	٧	٧	٧.	٧,	levels of exudate. Wicking spreads the fluid		ELA823
Form dressing	Mepilex Border- 1st choice	15 x 17.5cm	£3.80	√ √	√ √	V	√ √	V	√ √	throughout allowing for effective exudate		ELA714
Form dressing	Mepilex Border- 1st choice Mepilex Border- 1st choice	17 x 20cm 10 x 20cm	£5.23 £2.96	V	٧	V V	V V	V V	V V	management.		ELA078 ELA713
Form dressing	Mepilex Border- 1st choice	10 X 20CIII	12.90	V	V	V	V	_ <u> </u>	V			ELA/15
Form dressing	Biatain Adhesive	Sacral 17 x 27cm	£2.32	٧	٧	٧	٧	٧	٧	Foam dressings can be used to absorb low to high levels of exudate. Wicking spreads the fluid	To be used on incontinent patients	ELA636
Form dressing	Biatain Adhesive	Sacral 23 x 23cm	£4.27	٧	٧	٧	٧	٧	٧	throughout allowing for effective exudate management.		ELA050
Form dressing	Biatain non- Adhesive	5 x 7cm	£1.33	٧		٧			٧	Foam dressings can be used to absorb low to high levels of exudate. Wicking spreads the fluid		ELA414
Form dressing	Biatain non- Adhesive	10 x 10cm	£2.30	٧		٧			٧	throughout allowing for effective exudate management.		ELA039
Hydrocolloid	Duoderm extra thin- 1 st Choice	5 x 10cm (10)	£0.71	٧	٧	٧	٧	٧	٧	A biodegradable dressing, non-breathable and self- adhering to the skin. Easy to apply and impervious to	For dry to lightly exuding wounds can be used as a	ELM317
Hydrocolloid	Duoderm extra thin- 1 st Choice	7.5 x 7.5cm (5)	£0.75	٧	٧	٧	٧	٧	٧	bacterium. Facilitate rehydration and autolytic debridement of dry, sloughy or necrotic wounds. Also		ELM311
Hydrocolloid	Duoderm extra thin- 1 st Choice	10 x 10cm(10)	£1.23	٧	٧	٧	٧	٧	٧	suitable promoting granulation tissue. These dressings are not recommended for heavily exuding or infected wounds. Hydrocolloid dressings should	protectant for peri-ulcer-tissue on NPWT	ELM050
Hydrocolloid	Duoderm extra thin- 1 st Choice	9x 15cm(10)	£1.68	٧	٧	٧	٧	٧	٧	extend at least 2cm beyond the edge of the wound. Not for use on diabetic foot ulcers	uld and.	ELM101
Hydrocolloid	Granuflex	10 x 10cm(10)	£2.49	٧		٧			٧	A biodegradable dressing, non-breathable and self- adhering to the skin. Easy to apply and impervious to bacterium. Facilitate rehydration and autolytic		ELM141
Hydrocolloid	Granuflex	15 x 15cm(10)	£4.73	٧		٧			٧	debridement of dry, sloughy or necrotic wounds. Also suitable promoting granulation tissue. These dressings are not recommended for heavily exuding	To be used on incontinent patients E For dry to lightly exuding wounds can be used as a protectant for peri-ulcer-tissue on NPWT E For dry to lightly exuding wounds. Not for use on diabetics E Not for use on heavily exuding wounds	ELM143
Hydrocolloid	Granuflex	20x 20cm	£8.49	٧		٧			٧	or infected wounds. Hydrocolloid dressings should extend at least 2cm beyond the edge of the wound. Not for use on diabetic foot ulcers		ELM145
Hydrogel	Cutimed gel	8g	£1.71	٧	٧	٧	٧	٧	٧	Is a clear, amorphous used in the treatment of sloughy and necrotic tissue (autolytic debridement). Not suitable for heavily exuding wounds as can lead to skin maceration/excoriation	Not for use on heavily exuding wounds	ELA615
Charcoal dressing	Clinisorb	10 x 10cm	£2.01	٧	٧	٧	٧	٧	٧	Dressings containing activated charcoal that is effective in absorbing chemicals released from malodorous wounds. These can be used as a primary or secondary dressings in the management of exuding, purulent and contaminated wounds (e.g.	Apply on a primary or a source to the description	ELV051
Charcoal dressing	Clinisorb	10 x 20cm	£2.67	٧	٧	٧	٧	٧	٧	fungating) complicated by bacterial infection and offensive odour. May adhere to wounds that have dried out in which case a non-adhesive contact layer will be required. Can be used as a secondary dressing over a primary dressing to suppress odour. Moisture inactivated charcoal.	Apply as a primary or secondary dressing	ELV053



						Availab	le to					
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
Island dressing	Softpore	6x 7cm	£0.04	٧	√	٧	٧	٧	٧			EIJ023
Island dressing	Softpore	10 x 10cm	£0.08	√	√	√	√ √	√ √	√ √	Used to exects a horrior against further traums		EIJ013
Island dressing Island dressing	Softpore Softpore	10 x 15cm 10 x 20cm	£0.11 £0.16	√ √	√ √	√	V V	V V	V V	Used to create a barrier against further trauma, absorbing exudate and wound fluid	Mepore/Reliance in use in the hospital	EIJ014 EIJ024
	Softpore	10 x 25cm	£0.18	V	V	V	V	۷ ۷	V	absorbing extuate and wound had		EIJ025
	Softpore	10 x 30cm	£0.23	٧	٧	<u> </u>	٧	٧	٧			EIJ026
Island dressing	Clearpore	6x 7cm	£0.08	٧	٧	٧	٧	٧	٧			ELW714
Island dressing	Clearpore	10 x 10cm	£0.14	٧	√	٧	٧	٧	٧	Used to create a barrier against further trauma,		ELW716
Island dressing	Clearpore	10 x 20cm	£0.29	√	٧	√	٧	٧	٧	absorbing exudate and wound fluid	C-View being used in the hospital	ELW719
Island dressing	Clearpore	10 x 25cm	£0.34	٧	٧	٧.	٧	٧	٧	absorbing exadute and wound had		ELW720
Island dressing	Clearpore Opsite - for orthopaedic and	10 x 30cm	£0.45	٧	٧	٧	٧	٧	٧			ELW721
Post-operative dressing	sensitive skin	5x 6.5cm	£0.20	٧		٧	٧		٧			ELW052
Post-operative dressing		8.5x 9.5cm	£0.66	٧		٧	٧		٧			ELW051
Post-operative dressing		8x 10cm	£1.59	٧		٧	٧		٧			
Post-operative dressing		10x 12cm	£0.87	٧		٧	٧		٧	Used to absorb, retain blood and surgical exudate	-	ELW090
Post-operative dressing		10x 20cm	£1.15	٧		٧	٧		٧			ELW091
Post-operative dressing		10x 25cm	£1.15	٧		٧	٧		٧			ELW092
Post-operative dressing		10x 30cm	£1.89	٧		٧	٧		٧			ELW045
Post-operative dressing	Opsite visible- immediately post-operative only	8x 10cm	£1.59	٧								ELW306
Post-operative dressing		10x 15cm	£1.72	٧								ELW227
Post-operative dressing		10x 20cm	£2.07	٧						Used to absorb, retain blood and surgical exudate		ELW228
Post-operative dressing		10x 25cm	£2.28	٧								ELW229
Post-operative dressing		10x 30cm	£3.96	٧								ELW307
Film dressing	Hydrofilm	10x 12.5cm	£0.42	٧	٧	٧	٧	٧	٧	A thin polyurethane coated with an acrylic adhesive,		ELW243
Film dressing	Hydrofilm	10x 15cm	£0.52	٧	٧	٧	٧	٧	٧	they do not absorb wound exudate but allow for		ELW246
Film dressing	Hydrofilm	12x 25cm	£0.86	V	V	V	v	1	V	wound visualisation. For use as a secondary dressing	Currently using C-view in the hospital which is more	ELW241
-	-			•	·		•	· ·		over hydrogels or as a primary dressing to protect	expensive	
Film dressing	Hydrofilm	15x 20cm	£0.96	٧	٧	٧	٧	٧	٧	fragile skin. Do not use to cover foam dressings as it		ELW247
Film dressing	Hydrofilm	20x 30cm	£1.60	٧	٧	٧	٧	٧	٧	reduces MVTR.		ELW248
Barrier product	Medi Derma 5	cream 28g	£3.12	٧	٧	٧	٧	٧	٧		This product had previously been used in the	
Barrier product	Medi Derma 5	cream 92g	£3.78	٧	٧	٧	٧	٧	٧	Used to hydrate clinically dry skin allowing moisture	community with success. Product evaluations carried	
Barrier product	Medi Derma 5	Spray 30ml	£5.65	٧	٧	٧	٧	٧	٧	to be locked in	out in the hospital	
Barrier product	Medi Derma 5	No-sting aerosol		٧	√	٧	٧	٧	٧			
Antimicrobial: Manuka	Activon	25g	£2.09	٧	٧	٧	٧	٧	٧		Suitable for cavities	ELZ069
Antimicrobial: Manuka	Actilite	10x 10cm	£1.00	٧	v	√	٧	٧	٧	Suitable for critically colonized or locally infected wounds. Systemic antibiotics are indicated in cases of		EJE042
honey Antimicrobial: Manuka							-			overt wound infection where classical signs are evident. The usage of Manuka honey dressings should Reinforced alginate dressing for use on moderately		
honey	Algivon plus	5x 5cm	£2.07			√			٧		ELS549	
honey	Algivon plus	10x 10cm	£3.55			٧			٧		exuding wounds. Cut to the size of wound	ELS550
Sequestering	Cutimed Sorbact swab	4x 6cm	£1.93	٧	٧	٧	٧	٧	٧	Using a chemical agent which utilizes hydrophobic interactions within 30 seconds microbes are	Do not use in combination with ointments/creams as	ELY212
Antimicrobial: Sequestering	Cutimed Sorbact swab	7x 9cm	£3.13	٧	٧	٧	٧	٧	٧	irreversibly bound to the dressing thus removing bacteria from the wound.	the binding effect is impaired	ELY213

						Availab	le to					
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
Antimicrobial lodine	Inadine	5x5cm	£0.43	٧	٧	٧	٧	٧	٧	Suitable for critically colonised or locally infected wounds. Systemic antibiotics are indicated in cases of	For lightly exuding wounds. Do not use if patient has known to iodine, is pregnant or breastfeeding. Use	EKB501
Antimicrobial lodine	Inadine	10 x 10cm	£0.63	٧	٧	٧	٧	٧	٧	overt wound infection where classical signs are evident.	under medical supervision in infants up to 6 months old and patients with thyroid disease.	
Antimicrobial lodine	lodoflex	5g (other sizes to be recommended by TVN)	£4.62	٧	٧	٧	٧	٧	٧	Suitable for critically colonised or locally infected wounds. Systemic antibiotics are indicated in cases of overt wound infection where classical signs are evident.	Not for sure on dry necrotic tissue or on patients with known sensitivity to any of its ingredients. So not use on children, pregnant or lactating women or people with thyroid disorders or rental impairment.150g of iodoflex can be applied per patient per week. A single application should not exceed 50g.	EKBOO7
Antimicrobial Solution	Prontosan	350ml	£4.53	٧	٧	٧	٧	٧	٧	Used to clean wounds and prevent infection	Wound irrigation solution which aids the removal of	ELY248
Antimicrobial Solution	Prontosan	40ml	£0.66	٧	٧	٧	٧	٧	٧	osea to clean frounds and prevent infection	biofilm and wound debris	
Pressure reducing pad	Kerrapro Dermal Gel Pad*	Strip 50 x2.5cm	£5.95	٧		٧			٧	Polymer gel-shaped pad that redistributes pressure	Do not use on broken skin. KerraPro can be reused on	FES9915
Pressure reducing pad	Kerrapro Dermal Gel Pad*	Strip 30 x5 x 0.3cm	£6.38	٧		٧			٧	while protecting and padding bony prominences.	the same individual. Wash with soap and water and pat	FES9914
Pressure reducing pad	Kerrapro Dermal Gel Pad*	Pad 10 x10x 0.3cm	£4.24	٧		٧			٧	Used to reduce the risk of pressure ulcers	dry before applying	FES9912
Sterile forceps	Blue forceps-Farla Medical (FGR373)			٧	٧	٧	٧	٧	٧			FGR373
Disposable sterile scissors				٧	٧	٧	٧	٧	٧			
Bandages	Hosicrepe	6cmx 4m	£0.47	٧	٧	٧	٧	٧	٧			ECA811
Bandages	Hosicrepe	7.5cmx 4.5m	£0.62	٧	٧	٧	٧	٧	٧	1		ECA030
Bandages	Hosicrepe	10cmx 4.5m	£0.76	٧	٧	٧	٧	٧	٧	7	Cotton Crepe light support bandage	ECA009
Bandages	Hosicrepe	15cmx 4.5m	£1.08	٧	٧	٧	٧	٧	٧	7		ECA059
Bandages	Peha-haft(Latex free)	5cmx 4.5m	£0.45	٧	٧	٧	٧	٧	٧			EBA133
Bandages	Peha-haft(Latex free)	8cmx 4m	£0.59	V	V	V	V	V	V	7		EBA134
Bandages	Peha-haft(Latex free)	10cmx 4m	£0.69	V	V	٧	1	V	v	7		EBA130
Bandages	Actico	8cmx 6m	£3.82	V	· ·	V V	1	<u> </u>	V V			EBA032
				V	 	V V	V V	<u> </u>	V V	┪		EBA016
Bandages	Actico	10cmx 6m	£3.96	•	<u> </u>	•	٧		•	4		
Bandages	Actico	12cmx 6m	£5.05	٧		٧	٧	<u> </u>	√			EBA033
Bandages	K-Soft	10cmx 3.5m	£0.56	√	٧	٧	٧	٧	٧	_	Sub-compression bandage wadding	EPA028
Bandages	K-Lite	10cmx 4.5m	£0.76	٧	٧	٧			٧		second layer of four bandage system	ECA100
Bandages	K-Plus	10cmx 8.7m	£2.03	٧	٧	٧			٧		Type 3s light compression bandage Third layer in K-four multilayer system	ECA162
Bandages	K-Flex	10cmx 6m	£3.48	٧	٧	٧			٧		Type 3a cohesive compression bandage	ECD018
Bandages	K-Four kits	18-25cm	£7.39	٧	٧	٧		<u> </u>	٧	_		EVN008
Bandages	K-Four kits	25-30cm	£7.75	٧	٧	٧			٧	_	K-four multilayer kits	EVN155
Bandages	K-Four kits	Above 30cm	£10.79	٧	٧	٧			٧			EVN156
Bandages	K-two (latex free)	18-25cm	£9.41	٧	٧	٧			٧		V turo/2 lavar system)	ECA236
Bandages	K-two (latex free)	25-32cm	£10.28	٧	٧	٧			٧	7	K-two(2 layer system)	ECA236
Bandages	K-two Reduced	18-25cm	£9.20	٧					٧		V h/2 January	ECA205
Bandages	K-two Reduced	25-32cm	£10.04	٧					٧		K-two(2 layer system)-reduced compression	ECA206
Bandages	Mollelast	4x 4cm	£0.34	٧	٧	٧			٧		Toe bandaging	EBA064
Bandages	Viscopaste		£4.86	٧	٧	٧	٧	٧	٧		Zinc paste bandage	EFA011
Tape	Hypafix	5cmx 5m	2 1.00	V	V	٧	V	v	V		Zine paste sandage	EHR033
	Hypafix	10cmx 10m		V	٧	, ,	V	V	V	1		EHR113
Tape	Hypafix	15cmx 10m		V V	V V	v 1	- 1	- 1	V V	1		EHR031
Tape	Hypafix	+		V V	V	V √	V V	V V	V	┪		EHR117
Tape		20cmx 10m		•	•	•	V	V	<u> </u>			CUVII/
Tape	Micropore	2.5cmx 5m		٧	V	√	V	V .	٧	4		
Tape	Micropore	5cmx 5m		V	V	٧	V	V	٧			



						Availab	le to					
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
					20 0							
Stockinette	Multipurpose											
Stockinette	Comfinette	Size 01 20m roll	£1.10	٧	٧	٧	٧	٧	٧			EGJ041
Stockinette	Comfinette	Size 12 20m roll	£1.79	٧	٧	٧	٧	٧	٧			EGJ042
Stockinette	Comfinette	Size 56 20m roll	£2.85	٧	٧	٧	٧	٧	٧			EGJ043
Stockinette	Comfinette	Size 78 20m roll	£2.65	٧	٧	٧	٧	٧	٧			EGJ044
Stackinatta	Comfinette	Size T1 LARGE LEGS	£4.84	٧	V	V	٧	-1	٧			EGJ047
Stockinette	Comfinette	20M ROLL	14.04	V	V	v	V	V	٧			EGJ047
Stockinette	Comfifast	Green line 10m roll	£3.90	٧	V	٧	V	V	V			
				V		1	V		V			50007
Stockinette	Comfifast	Blue line 10m roll	£3.31	V	V	v	V	V	V			EGP007
Stockinette	Comfifast	Yellow line 10m roll	£4.85	٧	٧	٧	٧	٧	٧			EGP008
Dressing packs	Softdrape	Small	£0.43	٧	٧	٧	٧	٧	٧			
Saline	Irripod saline		0.22/each	٧	٧	٧	٧	٧	٧			
Skin closure	Steristrips	6x 75mm3strips	£11.67	٧	٧	٧	٧	٧	٧			EIR503
Skin closure	Steristrips	3x 75mm5strips	£11.87	٧	٧	٧	٧	٧	٧			EIR501
Probe	Probe-for cavity wounds		£0.45	٧	V	٧	٧	٧	٧			
	Kerramax - Highly Exuding											
Superabsorbent dressing	wounds	10 x 10cm	£1.37						٧			EME023
	Kerramax - Highly Exuding											
Superabsorbent dressing	wounds	10 x 22cm	£1.80						٧			EME024
	Kerramax - Highly Exuding											
Superabsorbent dressing	wounds	20 x 30cm	£3.62						٧			EME025
	Sorbion Sachet S - Highly											
Superabsorbent dressing		10 x 10cm	£1.49						٧			EJE156
	Exuding wounds Sorbion Sachet S - Highly											-
Superabsorbent dressing		20 x 10cm	£1.99						V			EJE157
	Exuding wounds											
Superabsorbent dressing	Sorbion Sachet S - Highly	20 x 20cm	£3.49						V			EJE158
	Exuding wounds											
Foam Dressing	Mepilex Lite	6 x 8.5cm	£1.90						٧		For Neonates or very fragile skin	ELA182
Foam Dressing	Mepilex Lite	10 x 10cm	£2.28						٧		, -	ELA184
Foam Dressing	Polymem	10 x 10cm	£2.87						٧		For use in radiotherapy reaction only	ELA303
Foam Dressing	Polymem	17 x 19cm	£8.32						٧		,	ELA306
Antimicrobial:	Atrauman AG	5 x 5cm	£0.59						v			EKB039
Silver dressing												
Antimicrobial:	Atrauman AG	10 x 10cm	£1.37						v		Must be removed prior to x-ray, MRI, Ultrasound Or	EKB040
Silver dressing	Acradinari	10 X 10cm	11.57						·		Diathermy	ENDOTO
Antimicrobial:	Atrauman AG	10 x 20cm	£2.37						v			EKB041
Silver dressing	Attauman AG	10 X 20CIII	12.57						·			LKDO41
Antimicrobial:	Atrauman AG + Extra	5 x 5cm	£1.81						v			
Silver dressing	Attautilati AG + Extra	J X JUII	11.61						٧		Absorbs wound fluid, disrupts and breaks down biofilm	
Antimicrobial:	Atrauman AG + Extra	10 x 10cm	£4.32						- 1		to expose and kill bacteria, while preventing biofilm	
Silver dressing	Attautilati AG + Extra	10 X 10CIII	14.52						٧		reformation.	
Antimicrobial:	Atrauman AG + Extra	15 x 15cm	£8.13						V		Should overlap 1cm onto skin surrounding the wound	
Silver dressing	Atrauman AG + Extra	15 X 15CM	18.13						v			
Antimicrobial:												
Silver dressing	Kerra Contact AG	5 x 5cm	£3.32						٧		An advanced silver wound contact layer dressing. The	ELY597
											dressing consists of three layers: two non-adherent	
Antimicrobial:	Kerra Contact AG	10 x 12.5cm	£8.11						v		polyethylene mesh wound contact layers and one	ELY598
Silver dressing	Nerra Contact / C	10 X 12.5011	20.11						·		absorbent core made of polyester. All three layers are	221330
Antimicrobial:											coated with Ag oxysalts	
Silver dressing	Kerra Contact AG	15 x 15cm	£14.22						V		coated with Ag oxysalts	ELY599
Antimicrobial: lodine	lodoflex	10g							V			EKB008
Antimicrobial: lodine	lodoflex	10g 17g							V V		Currently being used	EKB008
			£5.37	 			 	-	V		Currently being used	EKB009 EKB012
Antimicrobial: lodine	lodosorb	10g ointment	13.3/						v		Currently being used	LNDU12
Antimicrobial: Solution	Prontosan Wound Gel X		£14.30						٧		Gel containing betaine (a gentle effective surfactant that penetrates, disturbs and removes biofilm and	
Antimicrobial: Solution	Prontosan Wound Gel X		£29.15						٧		wound debris) and PHMB (to help control bacterial levels on the wound).	
			<u> </u>	<u> </u>				<u> </u>				



						Availab	le to					
Dressing Category	New Formulary Product	Sizes	Prices	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use:	Comment	Code
Antimicrobial: Solution	Flaminal hydro	15g	£8.11						٧		Alginate gel containing two antimicrobial enzymes, glucose oxidase and lactoperoxidase, which exert their effect without damaging healthy skin cells. Debrides the wound and manages moisture balance. Should not be used where patients have a previous sensitivity reaction to alginate dressings or to polyethylene glycol.	ELGUZI
Antimicrobial: Solution	Flaminal forte	15g	£8.10						٧		Same contra-indications as above. Also should not be	ELG022
Antimicrobial: Solution	Flaminal forte	50g	£27.61						٧		used in full-thickness (third-degree)burns	ELG023
Barrier product	Proshield	Barrier cream 115g	£11.73						٧		A dimethicone-based skin protectant containing	ELG023
Barrier product	Proshield	Spray cleanser 235ml	£7.44						٧		copolymer bio-adhesives. Indicated for use on both intact skin and skin that is damaged as a result of incontinence.	
Larvae	Biomonde								V		For Hospital use ONLY Must not be used on wounds that have a tendency to bleed or are close to a major blood vessel. Should not be used on patients on anticoagulants where the relevant clotting maker is not within clinical range. Should be used with caution on wounds over/adjacent to exposed organs or leading to a body cavity and then only under close supervision.	
Protease modulating ointment	Cadesorb	10g	£5.88						٧		Protease modulating ointment that reduces local	ELZ085
Protease modulating ointment	Cadesorb	20g	£5.98						٧		wound pH to around 5 and thereby modulates protease activity.	ELZ319
Debriding pad	Debrisoft	10x 10cm	£7.11						٧		For removal of superficial slough and debris including biofilm from wounds and skin during wound bed preparation.	ELZ354
Hydrogel	KerraLite Cool	6x 6cm	£1.99						٧		Pro-ionic hydrogel contact layer Contra-indicated for	
Hydrogel	KerraLite Cool	8x 8cm	£2.30						٧		thickness wounds; or as a covering for deep, narrow	
Hydrogel	KerraLite Cool	11x 11cm	£3.07						٧		cavities or sinuses	
Negative pressure consumables	Renasys	Sml foam dressing kit	£25.38						٧			
Negative pressure consumables	Renasys	Med foam dressing kit	£28.23						٧			
Negative pressure consumables	Renasys	Lrg foam dressing kit	£33.42						٧			
Negative pressure consumables	Renasys	300ml canister	£21.26						٧			
Negative pressure consumables	Renasys	Go carry bag	£16.70						٧			
Negative pressure	Renasys	Go carry strap	£5.57						٧			
consumables Negative pressure	Renasys	Y connector with	£3.25						٧			
consumables Negative pressure	Renasys	30cm tubing Soft port Y connector	£3.58						٧			
consumables Negative pressure consumables	Renasys	7cmx 10cm hydro adhesive gel dressing	£34.64						٧			
Negative pressure consumables	PICO*	single use,disposable,no canister kit.2 dressings,1pump & 2batteries							٧			



Dressing Category	New Formulary Product	Sizes	Prices			Availab	le to			Suitable for use:	Comment	Code
Diessing Category	New Formulary Product	31263	Files	Hospital	LU Clinics	DN	PN	NH	TVN	Suitable for use.	Comment	Code
Negative pressure consumables	PICO*	10x 20cm	£137.87						٧			ELZ653
Negative pressure consumables	PICO*	10x 30cm	£137.87						٧			ELZ654
Negative pressure consumables	PICO*	10x 40cm	£139.45						٧			ELZ650
Negative pressure consumables	PICO*	25x 25cm	£139.85						٧			ELZ481
De-sloughing dressing	Urgoclean	6 x6	£0.73						٧		A sterile, non-adherent slough trapping poly-absorbent fibre dressing with the TLC healing matrix to promote wound healing and enable pain-free dressing changes.	ELZ404
De-sloughing dressing	Urgoclean	10 x10 cm	£1.77						٧			ELZ405