



Medicine Supply Notification

MSN/2024/032

Insulatard[®] InnoLet[®] (insulin isophane human) 100units/ml suspension for injection
3ml pre-filled disposable devices

Tier 2 – medium impact*

Date of issue: 21/03/2024

Link: [Medicines Supply Tool](#)

Summary

- Insulatard[®] InnoLet[®] (insulin isophane human) 100units/ml suspension for injection 3ml pre-filled disposable devices are being discontinued with remaining stock exhausted by the end of May 2024.
- Humulin[®] I KwikPen[®] (insulin isophane human) 100units/ml suspension for injection 3ml pre-filled pens remain available and can support increased demand.
- Insulatard[®] Penfill[®] and Insulatard[®] vials (insulin isophane human 100units/ml) remain available but cannot support any increased demand.

Actions Required

Prescribers should:

- not initiate patients on Insulatard[®] InnoLet[®]
- identify and review patients currently prescribed Insulatard[®] InnoLet[®];
- consider prescribing Humulin[®] I KwikPen[®] pre-filled pens where appropriate, taking into account the patient's manual dexterity, vision, ability to use new device correctly, and whether support is required with administering the dose; and
- seek advice from specialist diabetes team on management options if the above option is not considered suitable.

Prescribers and pharmacies should:

- ensure that all patients initiated on a new device are counselled on the change in device, provided with appropriate needles, training on use of new device, including signposting to training videos, and advised closer monitoring of blood glucose levels may be required (see Supporting Information).

*Classification of Tiers can be found at the following link:

<https://www.england.nhs.uk/publication/a-guide-to-managing-medicines-supply-and-shortages/>

Supporting information

Human isophane insulin

Onset of action is within 1½ hours, reaches a maximum effect within 4–12 hours and the entire duration of action is approximately 24 hours.

InnoLet® delivery device

This was developed for people with a visual or physical impairment to be able to set and inject insulin doses. It is a compact pre-filled pen able to deliver 1 to 50 units in increments of 1 unit. It comprises a pre-filled insulin reservoir and a large simple to use dial that has larger print numbers and gives an audible click each time the dial is moved. The dial is turned to the number of units to be delivered and the top of the device is pressed to administer the dose, after which the dial returns to zero.

Humulin® I KwikPen®

A pre-filled / disposable pen injector containing a 3ml cartridge which delivers up to 60 units per dose in single unit increments.

Links to further information

[SmPC: Insulatard® InnoLet® pens](#)

[Package leaflet: Insulatard® InnoLet®](#)

[SmPC: Humulin® I KwikPen® pens](#)

[BNF: Isophane insulin](#)

[BNF: Insulin preparations](#)

[NICE guidance: Type 1 diabetes in adults \(insulin therapy\)](#)

Enquiries

If you have any queries, please contact DHSCmedicinesupplyteam@dhsc.gov.uk.